

Meta.VOC

Proprietary bonded and crosslinked phase.


- Developed for analysis of volatile organic compounds (VOC)
- Intermediate polarity column

Meta.VOC Equivalent Phase

Agilent: DB-502.2, HP-VOC
Supelco: VOCOL
Restek: Rtx-502.2

Meta.VOC

Internal Diam.(mm)	Length (m)	Film Thickness (µm)	Temp limits (°C)	Part. N° (P/N)
0,18	20	1,00	-20 to 240/250	TR-941084
0,20	10	1,20	-20 to 240/250	TR-941249
0,25	30	1,50	-20 to 240/250	TR-941532
	60	1,50	-20 to 240/250	TR-941562
0,32	60	1,80	-20 to 240/250	TR-941863
	60	3,00	-20 to 230/240	TR-943063
0,53	30	3,00	-20 to 230/240	TR-943035
	60	3,00	-20 to 230/240	TR-943065
	105	3,00	-20 to 230/240	TR-9430K5



TRB-608

Proprietary bonded and crosslinked phase.

- Specifically designed for analysing chlorinated pesticides and PCBs
- Designed for the EPA 508, 608 and 8080 methods.

TRB-608 Equivalent Phase

Agilent: HP-608
Supelco: SPB-608
SGE: BP-608

TRB-608

Internal Diam.(mm)	Length (m)	Film Thickness (µm)	Temp limits (°C)	Part. N° (P/N)
0,18	20	0,18	-20 to 300/310	TR-360984
0,25	30	0,25	-20 to 300/310	TR-360232
0,53	15	0,50	-20 to 290/300	TR-360515
	30	0,50	-20 to 290/300	TR-360535

TRB 608

Column: **TRB-608**, P/N TR-360232
Dimensions: 30m x 0.25mm x 0.25mm
Injection: 1 µL Grob Test, split 1:100, 280°C
Carrier gas: H2, constant pressure 11 psi (75.8 Kpa)
Oven temperature: 40°C @ 6°C/min to 300°C(5min)
Detector: FID, 340°C

