

Teknokroma Headspace Vials



Headspace Vials



Cat.No	Description	Pk
TR-400085	10 ml Clear Glass, Beveled Edge, Flat bottom, Crimp HS Vial	100
TR-400086	10 ml Clear Glass, Beveled Edge, Rounded-Flat bottom, Crimp HS Vial	100
TR-400037	20 ml Clear Glass, Beveled Edge, Flat bottom, Crimp HS Vial	100
TR-400087	20 ml Clear Glass, Beveled Edge, Rounded-Flat bottom, Crimp HS Vial	100
TR-400088	20 ml Amber Glass, Beveled Edge, Rounded-Flat bottom Crimp HS Vial short neck (PE)	100



Cat.No	Description	Pk
TR-400030	20 mm Open Top Aluminium Crimp Cap (10 mm Hole) PTFE/Natural Silicone Septa	100
TR-400038	20 mm Open Top Magnetic Crimp Cap (8 mm Hole) Red & Silver PTFE/Silicone Septa	100
TR-400059	20 mm Open Top Magnetic Crimp Cap (8 mm Hole) Blue & Silver PTFE/Silicone Septa	100

20 mm Crimp Seals with Prefitted Septa for Headspace Vials



Cat.No	Seal Type	Septa	Pk
CC-C4020-34A	Regular	Ivory PTFE/Red Rubber, (-40 to +100°C)	100
CC-C4020-34AP	Pressure Release	Ivory PTFE/Red Rubber, (-40 to +100°C)	100
CC-C4020-39A	Regular	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-43A	Magnetic Steel	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-43AP	Pressure Release	Clear PTFE/Gray Butyl (-40 to +120°C)	100
CC-C4020-36A	Regular	Gray PTFE/Black Molded Butyl (-40 to +125°C)	100
CC-C4020-36AP	Pressure Release	Gray PTFE/Black Molded Butyl (-40 to +125°C)	100
CC-C4020-32A	Regular	Tan PTFE/White Silicone (-60 to +200°C)	100
CC-C4020-32AP	Pressure Release	Tan PTFE/White Silicone (-60 to +200°C)	100
CC-C4020-42A	Magnetic Steel	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100
CC-C4020-42AP	Pressure Release	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100
CC-C4020-37AP	Pressure Release	Aluminium Foil/White High Temp. Silicone (-60 to +220°C)	100
CC-C4020-31	Molded Polypropylene Storage Cap	SepCap Integral Molded Polypropylene	250
CC-C4020-42A	Magnetic Steel	Clear PTFE/Translucent Blue Silicone (-60 to +200°C)	100

Crimpers, Dicers, and Decapping Pliers for 20 mm Top Vials

Cat.No	Description	Use	Pk
CC-C4020-100	Manual Crimper	Attaches 20 mm aluminium crimp seals	1
CC-C4020-101	Decapping Pliers	Removes 20 mm aluminium crimp seals	1
CC-C4020-102	Manual Decrimper	Removes 20 mm aluminium seals without vial damage	1

CC-C4020-102

